

HRVT training laptops preparation for decrease group

1.Factory Reset

- a. After all backup, you need to do factory reset to clean all personal information in it.
- b. Shift key+Restart or go to Setting-Update & Security-Recovery-Advanced startup
- c. Select-Troubleshoot-Reset this PC-Remove everything-Reset –(if asked,choose keep my computer, then it will take shorter time)
- d. It will take several hours (1~2 hours) to reset and install Windows
- e. When it initialize the computer, do not connect network, just leave it offline.
- f. Make account name as “HRVT_#”

2.Removing all unnecessary functions and icons

- a. Unpin “Store”, “Amazon”, or any other icons in Taskbar
- b. Uninstall McAfee : Settings-System-App-uninstall completely
- c. Turn off Notification: Settings-System-Notifications and actions-turn off all
- d. Change sleep time: Settings-System-Power & Sleep-select 1 hour for every option
- e. **Search “control panel” in the search box in left on bottom task bar. And go to “System and Security”-“Sytem”-“Advanced system setting”.**
- f. **Choose the tab “Hardware”- “Device Installtion Settings” and change Yes to No.**

3.App Installation

• Decrease group

- 1) Removing all unnecessary functions and icons
 - a. Unpin “Store”, “Amazon”, or any other icons in Taskbar
 - b. Uninstall McAfee : Settings-System-App-uninstall completely
 - c. Turn off Notification: Settings-System-Notifications and actions-turn off all
 - d. Change sleep time: Settings-System-Power & Sleep-select 1 hour for every option
 - e. **Search “control panel” in the search box in left on bottom task bar. And go to “System and Security”-“Sytem”-“Advanced system setting”.**
 - f. **Choose the tab “Hardware”- “Device Installtion Settings” and change Yes to No.**
- 2) Display the instruction at Desktop Background
 - a. Instruction1.png
 - b. Settings-Personalization-Background (Picture)-Choose your picture (Browse and find the instruction file)-Choose a fit (Fit)
- ~~3) Copy the Music folder into Desktop.~~
- 4) Copy the installer folder and past it in the Downloads (after installing, remove this folder).
 - a. Decrease group: cgroup installer
- 5) Decrease group
 - a. Install emwave pro with license code(from emWave_license.txt)

- i. When it asks the account, you can put the device name “HRVT_#” as local account.
 - ii. Double click the emWave Pro in Desktop
 - iii. Go to edit –preference-Advanced and Check the External Games Communication Settings and TCP Port number 20480, Check Send Pulse Data.
 - iv. Go to edit –preference-General-Session Length should be 60.
 - v. **Go to edit-preference-General-Application Updates-uncheck “check for updates at start up”**
- b. Install “emWave_Pilot_1.0.12.msi”
 - i. If there is a previous version, need to remove all folders
 - 1. Remove c:\Program Files (x86)\emWave_Pilot\emWave_Pilot
 - 2. Remove C:\Users\HRVT_#\AppData\Roaming\emWave_Pilot
 - ii. **Then install new version (1.0.12)**
 - iii. Target folder: c:\Program Files (x86)\emWave_Pilot\emWave_Pilot
 - iv. Go to this folder and run the emWave_pulse_PC_App.ext
 - v. Click Start button and it will show error, but it should be find. But this first launching makes the data folder and info.csv file under the folder
(C:\Users\HRVT_#\AppData\Roaming\emWave_Pilot\Info,
C:\Users\HRVT_#\AppData\Roaming\emWave_Pilot\Log)
 - vi. ~~Copy the “Transfer” icon to Desktop~~
- c. Powershell Allow script execution
 - i. Run powershell (click the Start icon, scroll to "Windows PowerShell", click the "Windows PowerShell" folder, Right-click the "Windows PowerShell" application-more-Run as administrator)
 - ii. In the PowerShell console that opens, type:
Set-Executionpolicy RemoteSigned -Scope CurrentUser
and hit the "Enter" key. You'll get a warning asking you if you're sure you want to do this. Type 'Y' and hit "Enter".
 - iii. Close the console
- d. Move Powershell script and make shortcut on Background
 - i. Copy the start-exp.ps1 file to the Documents folder on the laptop.
 - ii. Copy the DailyTraining icon on Desktop
- e. Final Desktop background has only **Musie, DailyTraining, Transfer, HRVsocial and instruction image**. Erase emWave Pro icon in Background if exist.
- f. Erase installer folder in Downloads folder
- g. Empty Recycle Bin on Desktop

6) change ID in info.csv file

- a. Go to C:\Users\Participant_#\AppData\Roaming\emWave_Pilot\Info
- b. make a shortcut (right click of the folder and send to Desktop-then move the icon in the lab_only folder)
- c. Change the user name in info.csv file in Info folder.
 - i. Right-click and choose Edit
 - ii. [Decrease group Change the **User as Participant ID#, time as 300 or any number, threshold as “8”** and save it. Time unit is second.